

## Eurofins Genoma NIPT performance data related to CE-IVD accreditation

### Background

For many years now, Non-Invasive Prenatal Testing (NIPT) has become part of clinical practice and is used as a screening test by thousands of pregnant women in order to avoid invasive tests such as amniocentesis and chorionic villus sampling as much as possible and reduce the risks associated with them. NIPT analysis has expanded over time, allowing the study of not only trisomies 21, 18 and 13 (common fetal trisomies) and sex chromosome aneuploidies (SCA), but also microdeletion/duplication syndromes and genome-wide analysis that includes rare autosomal aneuploidies (RAAs) and partial deletions/duplications (1,2). Although, rarer fetal anomalies and partial deletions/duplications have a low prevalence, these anomalies can cause serious pregnancy complications and affect patient care. These complications include early miscarriage, fetal demise, intrauterine growth restriction, and birth defects (4).

NIPT based on massively parallel whole-genome sequencing, using paired-end sequencing (4), analyses small fragments of placenta- derived cell free DNA that are circulating in a pregnant woman's blood, allowing the detection of common aneuploidies, RAAs and deletions/duplications with high accuracy.

VeriSeq NIPT Solution V2 is a CE-IVD integrated platform and software that uses a CE-IVD protocol for genome-wide screening providing information on partial deletions/duplications >7Mb for all the autosomes and the aneuploidy status for all chromosomes.

Until now the potential to improve the performance of VeriSeq V2 protocol remained to be investigated.

### Scope and design of the study

In this study, we aimed to validate the clinical utility of Eurofins Genoma Group's NIPT performances in the detection of genome-wide fetal anomalies including common trisomies, SCAs, RAAs, and partial deletion/duplication >7 Mb and <7 Mb.

We **retrospectively analysed 71,883** patients undergoing NIPT from the general pregnancy population collected over a 24-month period between December 2019 and December 2021. We applied **Illumina VeriSeq NIPT solution v2** in combination with our **proprietary algorithm, developed in-house (Fig.1)**




Fig.1 Analysis and interpretation of the results with Eurofins Genoma's algorithm

#### ROMA

**Laboratori e Studi Medici**  
Sede Legale e Laboratorio  
di Ricerca e Sviluppo  
in Genetica Molecolare  
Via Castel Giubileo, 11  
00138

**Laboratorio Genetica  
Medica e Diagnostica  
Molecolare**  
Prelievi e Consulenze  
Via Castel Giubileo, 62  
00138

#### MILANO

**Laboratorio Genetica  
Molecolare e Studi Medici**  
Via Enrico Cialdini, 16  
(Affori Centre)  
20161

#### FIRENZE

**Laboratorio e studi medici**  
Via Cavour, 168r  
50121



The results of NIPT were compared with those of the clinical outcome. Clinical outcomes, i.e., clinical truth, for study cases were determined by invasive prenatal diagnostic techniques (cytogenetic analysis after chorionic villus sampling (CVS) and/or amniocentesis), as well as by ultrasound and new-born physical exam. Positive NIPT results for fetal aneuploidy were considered confirmed when validated by either invasive prenatal diagnostics or an anomaly observed on ultrasound that matched the positive NIPT call.

## Results

The results of this study show that the already excellent performances validated by Illumina are further improved when combined with the **test algorithm** designed by **our bioinformaticians** and **the experience of our professionals** allowing us to achieve an overall **very high sensitivity (99.49%) and specificity (99.88%)** (Table 1).

Table 1. NIPT Performances for common aneuploidies, SCAs and other abnormalities in 71883 pregnancies

Positive cases n=1011 Total follow-ups n=868	Triomy 21	Trisomy 18	Trisomy 13	SCA	Other abnormalities*	Overall performances
True positives	437	93	37	156	58	781
False positives	3	1	8	17	54	83
True negatives	71392	71775	71828	65598	46577	70872
False negatives	2	0	0	1	1	4
<b>Sensitivity (95% CI)</b>	<b>99.54%</b> (98.36% - 99.94%)	<b>99.9%</b> (96.11% - 100.00%)	<b>99.9%</b> (90.51% - 100.00%)	<b>99.36%</b> (96.50% - 99.98%)	<b>98.31%</b> (90.91% - 99.96%)	<b>99.49%</b> (98.70% - 99.86%)
<b>Specificity (95% CI)</b>	<b>99.9%</b> (99.99% - 100.00%)	<b>99.9%</b> (99.99% - 100.00%)	<b>99.99%</b> (99.98% - 100.00%)	<b>99.97%</b> (99.96% - 99.99%)	<b>99.88%</b> (99.98% - 99.99%)	<b>99.88%</b> (99.86% - 99.91%)
<b>PPV (95% CI)</b>	<b>99.32%</b> (97.92% - 99.78%)	<b>98.94%</b> (92.91% - 99.85%)	<b>82.22%</b> (69.82% - 90.24%)	<b>90.17%</b> (85.08% - 93.66%)	<b>51.79%</b> (45.08% - 58.42%)	<b>90.39%</b> (88.36% - 92.11%)
<b>NPV (95% CI)</b>	<b>99.9%</b> (99.99% - 100.00%)	<b>99.9%</b> (99.99% - 100.00%)	<b>99.9%</b> (99.99% - 100.00%)	<b>99.9%</b> (99.99% - 100.00%)	<b>99.9%</b> (99.99% - 100.00%)	<b>99.99%</b> (99.99% - 100.00%)

CI: Confidence Intervals; SCA: Sex Chromosome Aneuploidies.

Positive cases without follow-up that have been excluded from the positives reported in Table 1 (n°): T21 (49); T18 (14); T13 (10); SCA (36); Other abnormalities (34).

\*Rare autosomal aneuploidies, segmental anomalies and microdeletions are included.

In particular, for **common aneuploidies** such as **trisomy 21, 18 and 13** overall **sensitivity and specificity** were **99.65%** and **99.99%**. **Sex chromosome aneuploidies** showed a high reliability for XXY (Klinefelter syndrome), XYY (Jacobs syndrome) and XXX (trisomy X) anomalies and a slightly lower reliability for monosomy X (Turner syndrome) (Table 2), confirming excellent overall **sensitivity (99.36%)** and **specificity (99.97%)** (Table 1).

### ROMA

Laboratori e Studi Medici  
Sede Legale e Laboratorio  
di Ricerca e Sviluppo  
in Genetica Molecolare  
Via Castel Giubileo, 11  
00138

Laboratorio Genetica  
Medica e Diagnostica  
Molecolare  
Prelevi e Consulenze  
Via Castel Giubileo, 62  
00138

### MILANO

Laboratorio Genetica  
Molecolare e Studi Medici  
Via Enrico Cialdini, 16  
(Affori Centre)  
20161

### FIRENZE

Laboratorio e studi medici  
Via Cavour, 168r  
50121



**Table 2. NIPT Performances for sex chromosome aneuploidies**

Sex chromosome aneuploidies	X0	XXX	XXY	XYY
True positives	52	27	51	26
False positives	13	0	3	1
True negatives	65724	65775	65747	65776
False negatives	1	0	0	0
Sensitivity (95% CI)	<b>98.11%</b> (89.93% - 99.95%)	<b>99.99%</b> (87.23% - 100.00%)	<b>99.99%</b> (93.02% - 100.00%)	<b>99.99%</b> (86.77% - 100.00%)
Specificity (95% CI)	<b>99.98%</b> (99.97% - 99.99%)	<b>99.99%</b> (99.99% - 100.00%)	<b>99.99%</b> (99.99% - 100.00%)	<b>99.99%</b> (99.99% - 100.00%)
PPV (95% CI)	<b>80%</b> (69.88% - 87.34%)	<b>99.99%</b> (99.99% - 100.00%)	<b>94.44%</b> (84.57% - 98.14%)	<b>96.3%</b> (78.55% - 99.46%)
NPV (95% CI)	<b>99.99%</b> (99.99% - 100.00%)	<b>99.99%</b> (99.99% - 100.00%)	<b>99.99%</b> (99.99% - 100.00%)	<b>99.99%</b> (99.99% - 100.00%)

CI: Confidence Intervals; SCA: Sex Chromosome Aneuploidies.

Positive cases without follow-up that have been excluded from the positives reported in Table 2 (n°): X0 (18); XXX (6); XXY (7); XYY (5).

The data obtained for rare trisomies reveal **high sensitivity and specificity (99.99%; 99.92%)**, although with a low positive predictive value due to the high rate of fetoplacental mosaicism and the risk of early spontaneous abortion found in these cases (Table 3). The clinical usefulness of detecting rare trisomies is confirmed in relation to the possible effects of fetoplacental mosaicism on fetal growth and on the occurrence of other pregnancy related complications, especially in the third trimester of gestation. It is also particularly useful in identifying uniparental disomy of chromosomes subject to imprinting, as in the case of trisomy 15 whose rescue has led to the identification of two cases of Prader-Willi syndrome.

**Our NIPT analysis shows very good performances for the testing of segmental anomalies with dimensions greater than 7 Mb (sensitivity 99.99%; specificity 99.97%)** and in this case too, the presence of false positives is attributable to fetoplacental mosaicisms, as well as to the presence of maternal benign neoplasms such as uterine fibroids (Table 3).

The analysis of such a large statistical sample also allowed us to show for the first time the performances of the microdeletion syndrome tests, with promising results considering the rarity of the conditions investigated (sensitivity 83.33%; specificity 99.99%) (Table 3).

**ROMA**

Laboratori e Studi Medici  
Sede Legale e Laboratorio  
di Ricerca e Sviluppo  
in Genetica Molecolare  
Via Castel Giubileo, 11  
00138

Laboratorio Genetica  
Medica e Diagnostica  
Molecolare  
Prelievi e Consulenze  
Via Castel Giubileo, 62  
00138

**MILANO**

Laboratorio Genetica  
Molecolare e Studi Medici  
Via Enrico Cialdini, 16  
(Affori Centre)  
20161

**FIRENZE**

Laboratorio e studi medici  
Via Cavour, 168r  
50121



Table 3. NIPT Performances for rare autosomal aneuploidies, segmental chromosomal abnormalities, and microdeletions

Other anomalies	RAA	Segmental anomalies (>7 Mb) <sup>°°§§</sup>	Microdeletions* (segmental anomalies <7 Mb) <sup>§§</sup>
True positives	33	20	5
False positives	36	16	2
True negatives	46630	46681	28743
False negatives	0	0	1
<b>Sensitivity (95%CI)</b>	<b>99.99%</b> (89.42% - 100.00%)	<b>99.99%</b> (83.16% - 100.00%)	<b>83.33%</b> (35.88% - 99.58%)
<b>Specificity (95%CI)</b>	<b>99.92%</b> (99.89% - 99.95%)	<b>99.97%</b> (99.96% - 99.99%)	<b>99.99%</b> (99.99% - 100.00%)
<b>PPV (95%CI)</b>	<b>47.83%</b> (39.81% - 55.96%)	<b>55.56%</b> (43.37% - 67.11%)	<b>71.43%</b> (37.40% - 91.27%)
<b>NPV (95%CI)</b>	<b>99.99%</b> (99.99% - 100.00%)	<b>99.99%</b> (99.99% - 100.00%)	<b>99.99%</b> (99.99% - 100.00%)

CI: Confidence Intervals; RAA: Rare Chromosomal Aneuploidies.

Positive cases without follow-up that have been excluded from the positives reported in Table 3 (n°): RAA (25); Segmental abnormalities >7Mb (7); Microdeletions (2).

\*\*Investigated microdeletions: Di George Syndrome, Cri-du-chat Syndrome, Prader-Willi Syndrome, Angelman Syndrome, 1p36 Deletion Syndrome, Wolf-Hirschhorn Syndrome, Jacobsen Syndrome, Langer-Giedion Syndrome, and Smith-Magenis Syndrome.

§§Details are showed in Suppl. Table 1 and 2

With regard to the occurrence of false-negative results, these were in line with the international scientific literature on the limits of NIPT. An in-depth analysis of the few cases tested (<0.005%) confirmed the essential need to combine the NIPT study with an accurate ultrasound evaluation for correct and complete prenatal assessment.

This study demonstrates that our NIPT-algorithm is reliable and accurate when applied to maternal DNA samples collected from pregnant women. The performance data, which are superior to those of an already highly reliable protocol such as Illumina VeriSeq protocol, confirm the crucial importance of **data interpretation**, guaranteed by the experience of our biologists, geneticists and bioinformaticians, and the importance of having access to a **large and constantly updated sample pool** for the design of a valid test algorithm.

#### ROMA

**Laboratori e Studi Medici**  
Sede Legale e Laboratorio  
di Ricerca e Sviluppo  
in Genetica Molecolare  
Via Castel Giubileo, 11  
00138

**Laboratorio Genetica  
Medica e Diagnostica  
Molecolare**  
Prelievi e Consulenze  
Via Castel Giubileo, 62  
00138

#### MILANO

**Laboratorio Genetica  
Molecolare e Studi Medici**  
Via Enrico Cialdini, 16  
(Affori Centre)  
20161

#### FIRENZE

**Laboratorio e studi medici**  
Via Cavour, 168r  
50121



**Bibliographic references:**

1. Fiorentino F, Bono S, Pizzuti F, Duca S, Polverari A, Faieta M, et al. The clinical utility of genome-wide noninvasive prenatal screening. *Prenat Diagn* 2017;37:593–601.
2. Ehrich M, Tynan J, Mazloom A, Almasri E, McCullough R, Boomer T, et al. Genome-wide cfDNA screening: clinical laboratory experience with the first 10,000 cases. *Genet Med* 2017;19:1332–7.
3. Bianchi, D. W., Platt, L. D., Goldberg, J. D., Abuhamad, A. Z., Sehnert, A. J., Rava, R. P., & Maternal Blood IS Source to Accurately diagnose fetal aneuploidy (MELISSA) Study Group (2012). Genome-wide fetal aneuploidy detection by maternal plasma DNA sequencing. *Obstetrics and gynecology*, 119(5), 890–901. <https://doi.org/10.1097/AOG.0b013e31824fb482>
4. Pertile, M. D., Flowers, N., Vavrek, D., Andrews, D., Kalista, T., Craig, A., Deciu, C., Duenwald, S., Meier, K., & Bhatt, S. (2021). Performance of a Paired-End Sequencing-Based Noninvasive Prenatal Screening Test in the Detection of Genome-Wide Fetal Chromosomal Anomalies. *Clinical chemistry*, 67(9), 1210–1219. <https://doi.org/10.1093/clinchem/hvab067>
5. Pertile, M. D., Halks-Miller, M., Flowers, N., Barbacioru, C., Kinnings, S. L., Vavrek, D., Seltzer, W. K., & Bianchi, D. W. (2017). Rare autosomal trisomies, revealed by maternal plasma DNA sequencing, suggest increased risk of fetoplacental disease. *Science translational medicine*, 9(405), eaan1240. <https://doi.org/10.1126/scitranslmed.aan1240>
6. Harasim, T., Neuhann, T., Behnecke, A., Stampfer, M., Holinski-Feder, E., & Abicht, A. (2022). Initial Clinical Experience with NIPT for Rare Autosomal Aneuploidies and Large Copy Number Variations. *Journal of clinical medicine*, 11(2), 372. <https://doi.org/10.3390/jcm11020372>
7. Dharajiya, N. G., Namba, A., Horiuchi, I., Miyai, S., Farkas, D. H., Almasri, E., Saldivar, J. S., Takagi, K., & Kamei, Y. (2015). Uterine leiomyoma confounding a noninvasive prenatal test result. *Prenatal diagnosis*, 35(10), 990–993. <https://doi.org/10.1002/pd.4629>
8. Klimczak, Amber M.a; Seli, Emrea,b; Scott, Richard T. Jra Noninvasive prenatal testing in women undergoing in vitro fertilization with preimplantation genetic testing, *Current Opinion in Obstetrics and Gynecology*: June 2021 - Volume 33 - Issue 3 - p 184-187 doi: 10.1097/GCO.0000000000000707
9. Hui, L., & Bianchi, D. W. (2020). Fetal fraction and noninvasive prenatal testing: What clinicians need to know. *Prenatal diagnosis*, 40(2), 155–163. <https://doi.org/10.1002/pd.5620>
10. Zhu, X., Lam, D., Chau, M., Xue, S., Dai, P., Zhao, G., Cao, Y., Cheung, S., Kwok, Y., Choy, K. W., Kong, X., & Leung, T. Y. (2020). Clinical Significance of Non-Invasive Prenatal Screening for Trisomy 7: Cohort Study and Literature Review. *Genes*, 12(1), 11. <https://doi.org/10.3390/genes12010011>
11. Heling, K. S., Chaoui, R., & Flöter, M. (2018). False-Negative NIPT Result for Trisomy 21. *Ultraschall in der Medizin (Stuttgart, Germany : 1980)*, 39(1), 92. <https://doi.org/10.1055/s-0043-120758>
12. Suzumori, N., Sekizawa, A., Takeda, E., Samura, O., Sasaki, A., Akaishi, R., Wada, S., Hamanoue, H., Hirahara, F., Sawai, H., Nakamura, H., Yamada, T., Miura, K., Masuzaki, H., Nakayama, S., Kamei, Y., Namba, A., Murotsuki, J., Yamaguchi, M., Tairaku, S., ... Sago, H. (2021). Retrospective details of false-positive and false-negative results in non-invasive prenatal testing for fetal trisomies 21, 18 and 13. *European journal of obstetrics, gynecology, and reproductive biology*, 256, 75–81. <https://doi.org/10.1016/j.ejogrb.2020.10.050>
13. Samura, O., & Okamoto, A. (2020). Causes of aberrant non-invasive prenatal testing for aneuploidy: A systematic review. *Taiwanese journal of obstetrics & gynecology*, 59(1), 16–20. <https://doi.org/10.1016/j.tjog.2019.11.003>

Roma, 14 July 2022

**ROMA**

Laboratori e Studi Medici  
Sede Legale e Laboratorio  
di Ricerca e Sviluppo  
in Genetica Molecolare  
Via Castel Giubileo, 11  
00138

Laboratorio Genetica  
Medica e Diagnostica  
Molecolare  
Prelievi e Consulenze  
Via Castel Giubileo, 62  
00138

**MILANO**

Laboratorio Genetica  
Molecolare e Studi Medici  
Via Enrico Cialdini, 16  
(Affori Centre)  
20161

**FIRENZE**

Laboratorio e studi medici  
Via Cavour, 168r  
50121

